

Hill-Rom® Transport Stretcher



Hill-Rom Transport Stretcher

Hill-Rom's transport stretcher (P8005) is a cost effective choice for Emergency, OR, PACU or ASC settings. This stretcher offers unsurpassed mobility, comfort and quality.

Mobility. Easy to steer and stop, this stretcher features an optional next generation Steering Plus™ system for crisp cornering and control. A standard redundant-locking, four-wheel braking system provided for a strong, stable platform.

Comfort. Patients remaining on stretchers for extended periods of time due to overcrowding are at risk for developing pressure ulcers. Our innovative, 5-inch air mattress—AccuMax Quantum™ Stretcher Surface—is the only non-powered, dynamic pressure redistribution stretcher mattress with Controlled Release™ Technology.

Quality. Our reinforced frames are made of heavy-gauge automotive-grade steel for long lasting quality. Hill-Rom is committed to helping caregivers manage an increasing obese population. Focused on patient and caregiver safety, Hill-Rom is the only company that designs and tests all stretchers to a weight capacity of 700 pounds.

At a glance

Hill-Rom Transport Stretcher key features:

- **700 lb weight capacity**
- **Low height** for safer stretcher ingress and egress
- **Therapeutic Support Systems:** Foam, Tempur-Pedic®, AccuMax Quantum*
- **Ergonomic push handles**
- **Integrated IV pole transport device** helps to maneuver stretcher and a portable IV pole with ease
- **Integrated O₂ holder and tray***
- **Redundant-locking, four-wheel brakes**
- **Power washable**
- **Three year parts warranty** and lifetime frame weld warranty

* Optional feature

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

Hill-Rom® Transport Stretcher

Features	Standard Features	Optional Features
700 lb weight capacity	●	
26" width	●	
30" width		●
Brake/Steer pedals – head and foot	●	
Brake/Steer pedals – all four sides		●
Corner steer caster	●	
Next generation Steering Plus™ system		●
Integrated O ₂ tank holder with utility tray		●
Ergonomic push handles		●
Ergonomic push handles with integrated IV pole transport holder		●
Permanent IV pole – 2 stage		●
Siderail system – fold down	●	
3" or 4" high density foam mattress with fire barrier	●	
5" Comfortline® multi-density foam with fire barrier		●
3" or 4" Tempur-Pedic® mattress		●
5" AccuMax Quantum™ stretcher pad (Elastimax or Hypolex)		●
Bumper colors – choice of red or light neutral bumpers with neutral bellows	Choice	

Accessories

- Accessory IV pole
- Ankle stirrups
- Armboard
- Chart holder
- Convertible footboard
- Detachable footboard
- Foot extender
- Infusion support system
- IV transporter
- Lateral x-ray cassette holder
- Liquid oxygen tank holder
- Mattress pad
- Oxygen tank holder
- Paper roll dispenser
- Patient tray
- Siderail covers
- Stretcher pillow
- Transport straps
- Utility tray

Technical Specifications

Hill-Rom® Transport Stretcher


Maximum weight limit.....	700 lbs (317.5 kg)
Low position – floor to top of surface.....	23" (58.4 cm)
High position – floor to top of surface.....	35" (88.9 cm)
Overall length.....	83" (210.8 cm)
Overall width.....	32" (81.3 cm)
.....or 36" (91.4 cm)	
Siderail length.....	47" (119.4 cm)
Siderail height above mattress/surface.....	14" (36.8 cm)
Mattress size.....	26" (66 cm) x 75" (190.5 cm)
.....or 30" (76.2 cm) x 75" (190.5 cm)	
Maximum head elevation.....	90°
Maximum Trend/Reverse Trend.....	18°
Floor to base clearance.....	3.5" (8.9 cm)
Caster size.....	8" (20.3 cm)



AccuMax Quantum™ VPC Mattress



Tempur-Pedic® Medical Mattress



5" Comfortline® Foam Mattress



O₂ Tank Holder



Grey push handles with IV Pole Transport Device



Blue ergonomic push handles with IV Pole Transport Device



Next Generation Plus Steering™ system provides crisp cornering and control



Integrated Utility Tray



Integrated 3 Stage IV Pole

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2011 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

145646 rev 4 2/8/11

USA 800-445-3730
Canada 800-267-2337

www.hill-rom.com