

Datex-Ohmeda Excel 210 Anesthesia System

**Heartland
Medical**
The O.R. Outfitters

Powerful Ventilation & Monitoring Capabilities

The Datex-Ohmeda Excel 210 Anesthesia system is a great machine for facilities where a small footprint is appreciated. The Excel 210 is an all-pneumatic alternative for close quarters.

- The unit has O₂, N₂O and medical Air standard
- Also standard are hypoxic guard & a GMS absorber
- Comes with optional 7000 or 7800 ventilator

About This Machine

» Physical Features

- Compact
- Differentiated pipeline cylinder gauges

» Operation

- Master On/Off switch; Activates electronics & pneumatics
- Link-25 hypoxic mixture guard system
- Pin-indexed, gas-specific dual O₂, Air and N₂O flow meters
- Selecatec vaporizer manifold for two Tech Series Vaporizers
- Guarded O₂ flush button
- Oxygen power outlet

» Properties

- Fixed shelving with added monitoring space
- Optional drawer storage cabinets
- Single-action brake/footrest
- Large, easy-running casters
- Full-height dovetail accessory mounting system



» Available with Optional Anesthesia Ventilator

- Accommodates 3 gases & 2 vaporizers



Heartland
Medical Sales & Services, LLC

www.heartlandmedical.com

info@heartlandmedical.com

1-502-671-1014 | 877-329-7775