



Heartland
Medical Sales & Services, LLC

Heartland LS 300 X
Multi-Port, Xenon Light Source

www.heartlandmedical.com



The Heartland LS 300 X is a brilliant, long-lasting and dependable light source with intensity available at the flip of a switch. With 1,000 hour lamp-life, specifically designed cooling fan and multi-port rotating turret, the LS 300 X brings the versatility and functionality you demand to the OR.

Lamp: 300-Watt Xenon
Color Temp: 6,000° Kelvin
Lamp Life: 1,000 Hours
Dimensions: 12.6" x 11.5" x 5.25" (320mm x 292mm x 133mm)
Weight: 11lbs / 5kg
Power: 100v - 240v 50/60 Hz
Compliance: CAN/CSA C22.2 No. 601.1
CAN/CSA C22.2 No. 601.1.2
UL 6061 -1
IEC/EN 60601-1
IEC/EN 60601-1-2

Heartland Medical Sales & Services, LLC
2601 Holloway Road Louisville, KY 40299 | Phone: 502-671-1014 | Fax: 502-671-1019
www.heartlandmedical.com