



Heartland
Medical Sales & Services, LLC

www.heartlandmedical.com

Heartland PS 8100 LED
10-Watt Portable LED Headlight System



The Heartland PS 8100 LED 10-Watt Portable LED headlight system was built for the operating room and delivers full functionality, versatility and intensity into the field. With the integrated brilliance of a 10-Watt, 6000° Kelvin LED and a lightweight headband, the PS 8100 LED is illuminating and comfortable.

Lamp: 10-Watt LED

Color Temp: 6,000° Kelvin

Lamp Life: 25,000 Hours

Battery Life: 4.5 Hours on High

Charge Time: 5 Hours

Pinpoint Spot Adjustment: 2.5" - .5" Spot at 12" 64mm - 13mm Spot at 305mm
3.5" - .75" Spot at 16" 89mm - 19mm Spot at 406mm

Power: 100v - 240v

Heartland Medical Sales & Services, LLC
2601 Holloway Road Louisville, KY 40299 | Phone: 502-671-1014 | Fax: 502-671-1019
www.heartlandmedical.com