

Advanced Anesthesia Solutions

Although its technical features are sophisticated, you'll be surprised how simple Apollo is to use. Once staff know how to operate one device from the Dräger family, they can use them all – thanks to the wide range of operating features common to all Dräger Medical equipment and geared to worldwide industry standards.

- Helps maximize O.R. utilization
- Speeds O.R. turnover
- Supports clinical decision making through enhanced monitoring capability
- Supports optimized workflow through best-in-class ergonomics



About this Machine

» Features

- Efficient, ergonomic design saves time and optimizes workflow
- Easy to use, easy to move, easy to lock into place with a central foot brake

» Operation

- Advanced ventilation capability provides ICU quality ventilation in the O.R.
- Spontaneous breathing support
- Electrically powered piston ventilator requires no drive gas
- Enhanced monitoring capability provides user with the information they need to make clinical decisions

» Anesthesia Ventilation Modes

- Manual, spontaneous, Volume Mode, Pressure Mode
- Optional/Synchronization: Pressure Support (PS), Volume Mode Autoflow

