



Heartland
Medical Sales & Services, LLC

HeartlandMedical.com

**STERIS
HARMONY 585**

Fully Refurbished

The Harmony 585 LED Surgical Light is a visualization solution that provides exceptional surgical site illumination while its HD-capable suspension renders virtually every image signal in brilliant clarity. Harmony 585 LED lights are also known to reduce energy consumption by 3/4 and give you a cooler, more comfortable room with 34% less beam heat. The Harmony lights have no bulb-change requirements.

FEATURES:

- Bright light - up to 41% more illuminance than other LED Lights
- Pure, white light - 4,400K color temperature \pm 300K
- A large, bright spot - 160,000 lux max
- 84 overlapping beams of light virtually eliminate shadows from the surgeons head and hands
- Better tissue discrimination with a CRI of up to 96
- Deep, saturated reds with an R9 value of up to 98
- Accommodates virtually all HD signals - regardless of manufacturer or signal type (DVI, HD-SDI, RGBHV-etc)
- Grows with you - add light-heads, ambient light, monitors and cameras
- Add Harmony Connect Point for video integration without costly or involved renovations
- Convert general rooms to specialty suites quickly and easily
- Count on years of reliability due to the unique cable management system
- Color temperature: 4,400K
- Max Central Illuminance: 160,000 lux (14,870 fc)
- Life: 30,000 hours
- Average pattern sizes: 7-11" (18-28cm)



HEARTLAND MEDICAL SALES & SERVICES, LLC

2601 HOLLOWAY ROAD LOUISVILLE, KY 40299 | 1-502-671-1014 | 877-329-7775

www.HeartlandMedical.com